

CRP Materials Provided

Materials	25 Test	50 Test	100 Tests	1000 Tests	
CRP Latex Reagent	1.0 ml	2.0 ml	2 x 2.0 ml	20 x 2.0 ml	
Reactive Control	0.5 ml	1.0 ml	1.0 ml	10 x 1.0 ml	
Nonreactive Control	0.5 ml	1.0 ml	1.0 ml	10 x 1.0 ml	
Glycine-Saline Solution	50 ml	100 ml	100 ml	10 x 100 ml	
0.05 ml Disposable Stirrer Pipets	25	50	100	1000	
Disposable Test Cards (6-well)	5	9	17	170	

ADDITIONAL MATERIALS REQUIRED

- Timing device
- 13 x 75 mm test tubes
- Volumetric pipet to deliver 0.25 ml
- Mechanical rotator (optional)



ARLINGTON SCIENTIFIC, INC.

(800) 654-0146 · (801) 489-8911 · FAX (801) 489-5552
info@arlingtonscientific.com www.arlingtonscientific.com