

FACE SHIELD

PPE Face Shield for Healthcare Workers
by Arlington Scientific

FDA Establishment Registration Number 1641328

Shield Specs

- High-quality eye and face protection
- Single use, one-piece face shield
- Optically clear, distortion-free PET film
- Double-sided anti-fog coating
- Fits easily over eye glasses
- Soft stretchable headband
- This product is NOT made with natural rubber latex

Shield Sizes

- ME-20633-014 Full size shield (9" x 13"), **25 pack**
- ME-20633-015 Full size shield (9" x 13"), **100 case***

*4 boxes of 25

Directions for Use:

The face shield is designed as a nonsterile, single use product intended for use as a shield over the face to avoid liquid splashing or splattering. The face shield should be properly disposed of after use.

Conforms to OSHA 29CFR 1910.1030 Bloodborne Pathogens (3) Personal Protective Equipment (x) Face Shield Standards.

This product does not meet the requirements of ANSI Z87.1-2015 standards for occupational/educational eye and face protection since it is not impact resistant nor does the shield provide dust protection.

