

RF Materials Provided

Materials	25 Test	100 Test	1000 Tests		
RF Direct Latex Reagent	1.0 ml	4.0 ml	10 x 4.0 ml		
Reactive Control	0.5 ml	1.0 ml	10 x 1.0 ml		
Nonreactive Control	0.5 ml	1.0 ml	10 x 1.0 ml		
Glycine-Saline Solution	50 ml	100 ml	10 x 100 ml		
0.05 ml Disposable Stirrer Pipets	25	100	1000		
Disposable Test Cards (6-well)	5	17	170		

ADDITIONAL MATERIALS REQUIRED

- Timing device
- Volumetric pipet to deliver 0.25 ml
- Mechanical rotator (optional)



ARLINGTON SCIENTIFIC, INC.

(800) 654-0146 · (801) 489-8911 · FAX (801) 489-5552
info@arlingtonscientific.com www.arlingtonscientific.com